

Product Datasheet

CYP3A4 peptide (orb217185)

Description	The peptide is used to block anti-Cytochrome P450 3A4/5 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CYP3A4; CYP3A5
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213829.
Purity	> 85%
Source	Synthetic
Uniprot ID	P20815, P08684
Entrez	1577, 1576
Expiration Date	12 months from date of receipt.