

Product Datasheet

CYP4A11 peptide (orb217186)

Description	The peptide is used to block anti-Cytochrome P450 4A11/22 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CYP4A11; CYP4A22
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213830.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q5TCH4, Q02928
Entrez	284541, 1579
Expiration Date	12 months from date of receipt.