

Product Datasheet

CYP8B1 peptide (orb217187)

Description	The peptide is used to block anti-Cytochrome P450 8B1 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CYP8B1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213831.
Purity	> 85%
Source	Synthetic
Uniprot ID	O88962, Q9UNU6
Entrez	1582, 13124
Expiration Date	12 months from date of receipt.