

Product Datasheet

CYP17A1 peptide (orb217189)

Description	The peptide is used to block anti-Cytochrome P450 17A1 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CYP17A1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213833.
Purity	> 85%
Source	Synthetic
Uniprot ID	P11715, P27786, P05093
Entrez	13074, 1586, 25146
Expiration Date	12 months from date of receipt.