

Product Datasheet

DAB1 peptide (orb217192)

Description	The peptide is used to block anti-DAB1 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rat
Tested Applications	Blocking
Target	DAB1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213836.
Purity	> 85%
Source	Synthetic
Uniprot ID	O75553, P97318, Q8CJH2
Entrez	266729, 1600, 13131
Expiration Date	12 months from date of receipt.