

Product Datasheet

DDIT3 peptide (orb217198)

Description	DDIT3 peptide
Reactivity	Human, Mouse, Porcine, Primate, Rat
Tested Applications	Blocking
Target	DDIT3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213842.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q62857 , P35639 , P35638
Entrez	13198 , 1649
Expiration Date	12 months from date of receipt.