

Product Datasheet

DIAPH2 peptide (orb217206)

Description	The peptide is used to block anti-DIAPH2 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	DIAPH2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213850.
Purity	> 85%
Source	Synthetic
Uniprot ID	O60879, O70566
Entrez	54004, 1730
Expiration Date	12 months from date of receipt.