

Product Datasheet

DIO3 peptide (orb217207)

Description	The peptide is used to block anti-DIO3 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	DIO3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213851.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q91ZI8, P49897, P55073
Entrez	1735, 29475, 107585
Expiration Date	12 months from date of receipt.