

Product Datasheet

DLX5 peptide (orb217209)

Description	The peptide is used to block anti-DLX5 antibody reactivity.
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	DLX5
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213853.
Purity	> 85%
Source	Synthetic
Uniprot ID	P50575, P70396, P56178
Entrez	25431, 13395, 1749
Expiration Date	12 months from date of receipt.