

Product Datasheet

DNM1 peptide (orb217210)

Description	The peptide is used to block anti-Dynamin 1
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	DNM1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213854.
Purity	> 85%
Source	Synthetic
Uniprot ID	P39053, Q05193, P21575
Entrez	13429, 1759, 140694
Expiration Date	12 months from date of receipt.