

Product Datasheet

ATN1 peptide (orb217215)

Description	The peptide is used to block anti-Atrophin-1 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	ATN1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213859.
Purity	> 85%
Source	Synthetic
Uniprot ID	P54259, O35126, P54258
Entrez	1822, 29515, 13498
Expiration Date	12 months from date of receipt.