

Product Datasheet

RCAN1 peptide (orb217216)

Description	The peptide is used to block anti-DSCR1 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	RCAN1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213860.
Purity	> 85%
Source	Synthetic
Uniprot ID	P53805 , Q6IN33 , Q9JHG6
Entrez	54720 , 1827 , 266766
Expiration Date	12 months from date of receipt.