

Product Datasheet

DSG3 peptide (orb217217)

Description	DSG3 peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	DSG3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213861.
Purity	> 85%
Source	Synthetic
Uniprot ID	O35902, P32926
Entrez	1830
Expiration Date	12 months from date of receipt.