

Product Datasheet

E2F4 peptide (orb217223)

Description	The peptide is used to block anti-E2F4 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	E2F4; E2F5
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213867.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8R0K9 , Q15329 , Q16254 , Q62814 , Q61502
Entrez	1874 , 104394 , 13559 , 1875
Expiration Date	12 months from date of receipt.