

Product Datasheet

EFNA4 peptide (orb217231)

Description	The peptide is used to block anti-Ephrin A4 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rabbit, Rat
Tested Applications	Blocking
Target	EFNA4
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213875.
Purity	> 85%
Source	Synthetic
Uniprot ID	O08542, P52798
Entrez	13639, 1945
Expiration Date	12 months from date of receipt.