

Product Datasheet

EIF2S1 peptide (orb217238)

Description	The peptide is used to block anti-EIF2S1
Reactivity	Human, Mouse, Porcine, Rat, Zebrafish
Tested Applications	Blocking
Target	EIF2S1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213882.
Purity	> 85%
Source	Synthetic
Uniprot ID	P05198, P68101, Q6ZWX6
Entrez	1965, 54318, 13665
Expiration Date	12 months from date of receipt.