

Product Datasheet

EIF4G2 peptide (orb217242)

Description	The peptide is used to block anti-EIF4G2 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rabbit, Rat
Tested Applications	Blocking
Target	EIF4G2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213886.
Purity	> 85%
Source	Synthetic
Uniprot ID	P78344, Q62448
Entrez	1982, 13690
Expiration Date	12 months from date of receipt.