

Product Datasheet

ELAVL1 peptide (orb217244)

Description	The peptide is used to block anti-HUR antibody reactivity.
Reactivity	Human, Porcine, Sheep
Tested Applications	Blocking
Target	ELAVL1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213888.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q15717
Entrez	1994
Expiration Date	12 months from date of receipt.