

Product Datasheet

ELF5 peptide (orb217246)

Description	The peptide is used to block anti-ELF5 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	ELF5
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213890.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8VDK3, Q9UKW6
Entrez	13711, 2001
Expiration Date	12 months from date of receipt.