

Product Datasheet

ELK1 peptide (orb217248)

Description	The peptide is used to block anti-ELK1 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	ELK1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213892.
Purity	> 85%
Source	Synthetic
Uniprot ID	P19419 , P41969
Entrez	2002 , 13712
Expiration Date	12 months from date of receipt.