

Product Datasheet

EPHA5 peptide (orb217255)

Description	The peptide is used to block anti-EPHA5 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	EPHA5
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213899.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q60629 , P54757 , P54756
Entrez	13839 , 2044 , 79208
Expiration Date	12 months from date of receipt.