

Product Datasheet

NR2F6 peptide (orb217258)

Description	The peptide is used to block anti-NR2F6 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	NR2F6
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213902.
Purity	> 85%
Source	Synthetic
Uniprot ID	P10588, O09017, P43136
Entrez	245980, 2063, 13864
Expiration Date	12 months from date of receipt.