

Product Datasheet

ETS1 peptide (orb217263)

Description	The peptide is used to block anti-ETS1
Reactivity	Human, Mouse, Porcine, Rabbit, Rat
Tested Applications	Blocking
Target	ETS1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213907.
Purity	> 85%
Source	Synthetic
Uniprot ID	P14921, P41156, P27577
Entrez	23871, 2113, 24356
Expiration Date	12 months from date of receipt.