

Product Datasheet

F10 peptide (orb217271)

Description	The peptide is used to block anti-Factor X HC antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	F10
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213915.
Purity	> 85%
Source	Synthetic
Uniprot ID	O88947, Q63207, P00742
Entrez	29243, 2159, 14058
Expiration Date	12 months from date of receipt.