

Product Datasheet

FGFR1 peptide (orb217283)

Description	The peptide is used to block anti-FGFR1
Reactivity	Human, Mouse, Porcine, Rabbit, Rat, Zebrafish
Tested Applications	Blocking
Target	FGFR1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213927.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q04589, P16092, P11362
Entrez	2260, 14182, 79114
Expiration Date	12 months from date of receipt.