

Product Datasheet

FGFR1 peptide (orb217284)

Description	The peptide is used to block anti-FGFR1
Reactivity	Human, Mouse, Porcine, Primate, Rabbit, Rat
Tested Applications	Blocking
Target	FGFR1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213928.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q04589 , P16092 , P11362
Entrez	2260 , 14182 , 79114
Expiration Date	12 months from date of receipt.