

Product Datasheet

FHL2 peptide (orb217289)

Description	The peptide is used to block anti-FHL2 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	FHL2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213933.
Purity	> 85%
Source	Synthetic
Uniprot ID	070433, Q14192, O35115
Entrez	14200, 2274, 63839
Expiration Date	12 months from date of receipt.