

Product Datasheet

FOS peptide (orb217301)

Description	The peptide is used to block anti-c-FOS antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat, Sheep
Tested Applications	Blocking
Target	FOS
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213945.
Purity	> 85%
Source	Synthetic
Uniprot ID	P01100, P12841, P01101
Entrez	14281, 314322, 2353
Expiration Date	12 months from date of receipt.