

Product Datasheet

FPR1 peptide (orb217303)

Description	The peptide is used to block anti-FPR1 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	FPR1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213947.
Purity	> 85%
Source	Synthetic
Uniprot ID	P33766, P21462
Entrez	14293, 2357
Expiration Date	12 months from date of receipt.