

Product Datasheet

GABRA4 peptide (orb217317)

Description	The peptide is used to block anti-GABRA4 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	GABRA4
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213961.
Purity	> 85%
Source	Synthetic
Uniprot ID	P28471, Q9D6F4, P48169
Entrez	2557, 140675, 14397
Expiration Date	12 months from date of receipt.