

## Product Datasheet

### GAPDH peptide (orb217323)

|                            |                                                              |
|----------------------------|--------------------------------------------------------------|
| <b>Description</b>         | The peptide is used to block anti-GAPDH antibody reactivity. |
| <b>Reactivity</b>          | Human, Mouse, Porcine, Primate, Rabbit, Rat                  |
| <b>Tested Applications</b> | Blocking                                                     |
| <b>Target</b>              | GAPDH                                                        |
| <b>Form/Appearance</b>     | Lyophilized powder                                           |
| <b>Storage</b>             | Shipped at 4°C. Store at -20°C for one year.                 |
| <b>Note</b>                | For research use only                                        |
| <b>Application notes</b>   | This peptide was used as immunogen to produce orb213967.     |
| <b>Purity</b>              | > 85%                                                        |
| <b>Source</b>              | Synthetic                                                    |
| <b>Uniprot ID</b>          | <b>P16858, P04406, P04797</b>                                |
| <b>Entrez</b>              | <b>24383, 2597, 14433, 100042025</b>                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                              |