

Product Datasheet

GATA4 peptide (orb217325)

Description	The peptide is used to block anti-GATA4 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	GATA4
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213969.
Purity	> 85%
Source	Synthetic
Uniprot ID	P43694, Q08369, P46152
Entrez	2626, 54254, 14463
Expiration Date	12 months from date of receipt.