

Product Datasheet

GBP3 peptide (orb217326)

Description	The peptide is used to block anti-GBP3 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	GBP3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213970.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9H0R5
Entrez	2635
Expiration Date	12 months from date of receipt.