

Product Datasheet

GCGR peptide (orb217328)

Description	Glucagon-like peptide 1 (GLP-1) receptor (GLP1R)/glucagon receptor (GCGR) dual agonist peptide that ...
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	GCGR
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213972.
Purity	> 85%
Source	Synthetic
Uniprot ID	P47871, P30082
Entrez	24953, 2642
Expiration Date	12 months from date of receipt.