

Product Datasheet

GCK peptide (orb217329)

Description	The peptide is used to block anti-Glucokinase antibody reactivity.
Reactivity	Human, Mouse, Rat, Sheep
Tested Applications	Blocking
Target	GCK
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213973.
Purity	> 85%
Source	Synthetic
Uniprot ID	P52792, P35557, P17712
Entrez	24385, 2645, 103988
Expiration Date	12 months from date of receipt.