

Product Datasheet

GDF1 peptide (orb217330)

Description	The peptide is used to block anti-GDF1 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	GDF1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213974.
Purity	> 85%
Source	Synthetic
Uniprot ID	P20863, P27539
Entrez	2657, 14559
Expiration Date	12 months from date of receipt.