

Product Datasheet

GDI1 peptide (orb217332)

Description	The peptide is used to block anti-GDI1 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rabbit, Rat, Sheep, Zebrafish
Tested Applications	Blocking
Target	GDI1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213976.
Purity	> 85%
Source	Synthetic
Uniprot ID	P31150, P50398, P50396
Entrez	25183, 14567, 2664
Expiration Date	12 months from date of receipt.