

Product Datasheet

GHRHR peptide (orb217333)

Description	The peptide is used to block anti-GHRH Receptor antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	GHRHR
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213977.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q02644 , Q02643 , P32082
Entrez	2692 , 14602 , 25321
Expiration Date	12 months from date of receipt.