

Product Datasheet

TCTN1 Protein Lysate (orb2173333)

Description	TCTN1 Protein Lysate
Preservatives	Lyophilized
Form/Appearance	Lyophilized
Storage	The lysate can be shipped at ambient temperature. Upon receiving, store the sample at -2°C. Lysate samples are stable for 12 months from date of receipt when stored at -2°C. Avoid repeated freeze-thaw cycles after reconstitution. Lysate samples can be reconstituted with SDS Sample Buffer. Prior to SDS-PAGE fractionation, boil the lysate for 5 minutes
Tag	C-MYC/DDK
Note	For research use only
MW	61.9 kDa
Uniprot ID	Q2MV58
NCBI	NP_078825
Expiration Date	6 months from date of receipt.