

Product Datasheet

GJB1 peptide (orb217335)

Description	The peptide is used to block anti-Connexin 32 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rat
Tested Applications	Blocking
Target	GJB1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213979.
Purity	> 85%
Source	Synthetic
Uniprot ID	P08034, P08033, P28230
Entrez	2705, 29584, 14618
Expiration Date	12 months from date of receipt.