

Product Datasheet

GLRB peptide (orb217340)

Description	The peptide is used to block anti-GLRB antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	GLRB
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213984.
Purity	> 85%
Source	Synthetic
Uniprot ID	P48167, P20781, P48168
Entrez	2743, 25456, 14658
Expiration Date	12 months from date of receipt.