

Product Datasheet

GNAI1 peptide (orb217342)

Description	The peptide is used to block anti-GNAI1 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	GNAI1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213986.
Purity	> 85%
Source	Synthetic
Uniprot ID	P10824, P63096, B2RSH2
Entrez	25686, 2770, 14677
Expiration Date	12 months from date of receipt.