

Product Datasheet

GP1BB peptide (orb217345)

Description	The peptide is used to block anti-CD42c antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	GP1BB
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213989.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9JJM7, P56400, P13224
Entrez	2812, 116727, 14724
Expiration Date	12 months from date of receipt.