

Product Datasheet

GPR18 peptide (orb217351)

Description	The peptide is used to block anti-GPR18 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	GPR18
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213995.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q14330 , A1A5S3 , Q8K1Z6
Entrez	2841 , 110168 , 679957
Expiration Date	12 months from date of receipt.