

Product Datasheet

GPR18 peptide (orb217352)

Description	The peptide is used to block anti-GPR18 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	GPR18
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213996.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q14330, A1A5S3, Q8K1Z6
Entrez	2841, 110168, 679957
Expiration Date	12 months from date of receipt.