

Product Datasheet

GPR19 peptide (orb217353)

Description	The peptide is used to block anti-GPR19 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	GPR19
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213997.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q15760, Q61121, P70585
Entrez	312787, 2842, 14760
Expiration Date	12 months from date of receipt.