

Product Datasheet

GPR20 peptide (orb217354)

Description	The peptide is used to block anti-GPR20 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	GPR20
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb213998.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q99678, Q8BYC4
Entrez	2843, 239530
Expiration Date	12 months from date of receipt.