

Product Datasheet

GRIA1 peptide (orb217363)

Description	The peptide is used to block anti-GLUR1
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	GRIA1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214007.
Purity	> 85%
Source	Synthetic
Uniprot ID	P19490, P23818, P42261
Entrez	2890, 50592, 14799
Expiration Date	12 months from date of receipt.