

Product Datasheet

GRM2 peptide (orb217367)

Description	The peptide is used to block anti-mGLUR2 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	GRM2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214011.
Purity	> 85%
Source	Synthetic
Uniprot ID	P31422 , P31421 , Q9QYS2 , Q14832 , Q14BI2 , Q14416
Entrez	2912 , 108069 , 2913 , 24415 , 108068 , 24416
Expiration Date	12 months from date of receipt.