

Product Datasheet

GRM4 peptide (orb217368)

Description	The peptide is used to block anti-mGLUR4 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	GRM4
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214012.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q68EF4 , P31423 , Q14833
Entrez	2914 , 24417 , 268934
Expiration Date	12 months from date of receipt.